

MANUAL VICAT APPARATUS (BS/EN)

Standard: BS EN 196-3/480-2

Determination of Setting Time & Consistency of Cement

The vicat frame consists essentially of a metal stand with a sliding rod. An adjustable indicator moves over a graduated scale. The needle or plunger is attached to the bottom end of the rod to make up the test weight of 300 g. The frame is supplied with basic BS EN accessories.

Technical Specifications :

Model Number	CL 3012 X / 003
Standard	BS/EN
Dimension (mm)	215(W)x160(L)x323(H)
Approx. Weight	4.0 kg


Set Consists Of :

Model Number	Parts Description	Qty
CL 3012 X / 003 – P 001	Conical Plastic Mould - Ø70/80x40 (H) mm	1no.
CL 3012 X / 003 – P 002	Initial Needle Ø 1.13 mm ± 0.05	1no.
CL 3012 X / 003 – P 003	Final Needle Ø 1.13 mm ± 0.05	1no.
CL 3012 X / 003 – P 004	Consistency Plunger Ø10mm	1no.
CL 3012 X / 003 – P 005	Glass Plate 4 x 98 x 98 mm	2nos.


*1 Copy of Instruction Manual

Optional Accessories :

Model Number	Parts Description
CL 3012 X / 003 – A 001	Additional Weight 700g (BS/EN)
CL 3012 X / 003 – A 002	Conical Brass Mould (BS)
CL 7034 X / 005 A	Glass Thermometer -10°C to 50°C x 1°C Div.



700 g Weight



Glass Thermometer



Conical Brass Mould